

University of Moratuwa

Department of Electronic and Telecommunication Engineering



BM4152

Bio-signal Processing

Multi-channel EEG: Analyzing using EEGLAB

K. P. Tharuka

200641T


Contents

Part 1: Installing EEGLAB and importing data.....	3
Part 2: Data preparation, time and frequency domain analysis.....	4
Part 3: Using independent component analysis (ICA) of EEG data	16
blinks	16
<i>i_LR</i>	20
<i>i_UD</i>	24
teeth-blink.....	28
noise-ON-blink	32
noise-OFF	36

Part 1: Installing EEGLAB and importing data

1. EEGLAB opensource MATLAB toolbox is installed according to the given instructions.
2. Then the *xdf import* plugin is installed and then using that the *trial_2.xdf* file is imported and named as *raw_data*.

Information of the dataset is as follows.



The screenshot shows the EEGLAB v2024.2 software window. The menu bar includes File, Edit, Tools, Plot, Study, Datasets, and Help. The main window displays the properties of a dataset named '#1: raw_data'. The properties listed are:

Property	Value
Filename:	none
Channels per frame	16
Frames per epoch	46896
Epochs	1
Events	42
Sampling rate (Hz)	256
Epoch start (sec)	0.000
Epoch end (sec)	103.184
Reference	unknown
Channel locations	No (labels only)
ICA weights	No
Dataset size (Mb)	3.5

Part 2: Data preparation, time and frequency domain analysis

3. Channel labels

3.1. Channel location information is given using the provided *electrode_locations.locs* file.


3.2. Then it shows in dataset information GUI as follows

Channel locations


Yes

3.3.

Channel locations used during the practical are as follows.




Channel locations in the *electrode_locations.locs* file.



By comparing the electrode positions in both cases, it is clear that the channel locations are correct and well aligned. However, small deviations may have occurred when placing the electrodes on the subject due to the irregular shape of the skull, which is not perfectly spherical. This could cause slight variations in the exact placement.

4. Visualize data in time domain.




5. Re-referencing

Re-referencing in EEG signal processing refers to the process of selecting a new reference point (or electrode) to improve signal quality and reduce the common mode noise. EEG signals are typically recorded as the difference in voltage between an active electrode and a reference electrode. By re-referencing, the recorded signals are recalculated relative to a different reference electrode or set of electrodes.

Average referencing is one of the common referencing techniques where the average of all electrodes is used as the new reference.

Here we can clearly observe that the signals have less noise and the some electrical activities are clearly visible.

Referenced and then normalized plot of EEG data is given below.



6. Event renaming


Events are renamed according to the information in practical guidance.

Table 3 - Sequence of events

Event key	Duration (s)	Description
-	10	Relaxed and eyes open (default)
A	5	
B	5	
A	5	
B	5	
-	5	Relaxed and eyes open (default)
C		
C		
C		
D		Center-left-right-center
-	5	Relaxed and eyes open (default)
E		Center-up-down-center
-	5	Relaxed and eyes open (default)
F	2	
-	5	Relaxed and eyes open (default)
F	2	
-	5	Relaxed and eyes open (default)
F, C	5	Together
-	5	Relaxed and eyes open (default)
G	5	
-	10	Relaxed and eyes open (default)
H	5	Near the right temporal lobe
H, C	5	
I	5	
H	5	Near the left temporal lobe
H, C	5	
I	5	

Event key	Event description	Event short-name
A	Eyes closed	i_close
B	Eyes open	i_open
C	Blink (close and open)	blinks
D	Eyeball movement left-right (eyes open)	i_LR
E	Eyeball movement up-down (eyes open)	i_UD
F	Teeth clenching	teeth
G	Head movement (random)	Head
H	Noise source ON	noise_ON
I	Noise source OFF	noise_OFF


time domain plot with renamed events



7. Segmenting data related to events

Method 1: based on events

Eye blinking events are segmented: Set the start time limit to -0.5 from the event and end time limit to 1.2 from the event.




Segmented dataset is saved as *blinks_trials*.

This method has few issues.









- Here we can observe that **some time windows include data from events before or after the target event.**
- In event-related data collection, a specific time window is typically chosen to capture data before and after an event. However, the accuracy of these markers is subjective, meaning it can vary based on judgment. This makes it **difficult to find an exact time window that always includes the intended event.** As a result, while the first three trials capture the blink, the last two miss it.



Method 2: based on time

In this method since we manually select the time stamps looking at the recorded data we can isolate the desired events.



Events	Time Stamps	Dataset name
i-close	[10.5 17; 23.5 32]	i closed
i-open	[17 23.5; 32 36]	i open
blinks	[36 38; 38 40; 40 42]	blinks
i-LR	[50.5 57.5]	i LR
i-UD	[66 73]	i UD
teeth	[80.5 83.5; 88.5 93.5]	teeth
teeth + blinks	[96 100.5]	teeth_blink
head	[109 118]	head
noise-OFF	[152.5 159; 180 183]	noise OFF
noise-ON + blinks	[145 149; 173 176]	noise ON blink







8.

All the datasets after completing the lab.


** noise_ON dataset is not prepared as it is not given in the question 7.2.4. but only for the question 10.6 the dataset is temporarily prepared.



9. Manual rejection of artifacts




Before




after

10.



I_closed



The spectral power shows a clear **alpha peak (8-12 Hz)**, which is characteristic of a relaxed brain state when the eyes are closed. This indicates dominant alpha rhythms, a common feature during rest. The **scalp maps** display concentrated alpha activity in the **occipital lobes** in that frequency band , confirming the brain's visual disengagement and relaxation during the eyes-closed condition.





i_LR




The EEG powers spectrum density plot shows low-frequency activity below 5Hz in the frontal region, which is a clear sign of eye movements. This is supported by the presence of horizontal equivalent current dipoles (ECDs) in scalp topography at 10Hz, indicating horizontal eye movements. The peak at around 10 Hz in channel 13 shows the activity in the occipital region, which suggests that the visual cortex is being stimulated.

I_UD




In contrast to i_LR here we can observe vertical equivalent current dipoles (ECDs) in scalp topography at 6Hz and 10Hz.

Teeth



The plots show a significant increase in spectral power, particularly in the higher frequency range (above 20 Hz). This spike in high-frequency activity is a clear indication of muscle movements, specifically from jaw clenching. Such muscle activity dominates the EEG recordings and overshadows the brain's typical electrical activity. The scalp maps confirm that the strong electrical signals are not due to brain activity but rather to muscle contraction, which creates noise in the EEG readings.

noise_ON




The spectral plots show a steep drop in power at lower frequencies, followed by a relatively flat curve with some small fluctuations. Typically, we would see peaks in the 5-30 Hz range, which are linked to brain rhythms. However, here the noise has overpowered these expected peaks

Part 3: Using independent component analysis (ICA) of EEG data


11.

blinks

Normalized time domain ICA Components of the blinks-dataset






Spectrum and topography of independent components related to blinks



Classification of ICA components


ICA Component	Category – eye, brain, muscle, power line noise, channel noise, other	Reason
1	Eye	Scalp topographies suggest Equivalent Current Dipoles near eyes, the power concentrated at low frequencies (below 5 Hz) and clearly visible large peaks in the time domain plot.
2	Eye	Similar to 1
3	Eye	Similar to 1
4	Eye/ Muscle	Scalp topography suggests Equivalent Current Dipoles near eyes, but some peaks are visible above 20Hz.
5	Other	Scalp topography suggests Equivalent Current Dipoles but power spectrum and time domain plot don't give a strong evidence .
6	Eye	Similar to 1 but time domain plot doesn't show large peaks as 1.
7	Brain	Scalp topography looks dipolar, Power spectrum decreases as frequency increases ($1/f$) most importantly a peak around 10Hz.
8	Eye	Similar to 1 but time domain plot doesn't show large peaks as 1.
9	Muscle	Power concentrated in higher frequencies than lower frequencies and Scalp topography localizes near ears.
10	Brain	Similar to 7
11	Muscle	Power concentrated in higher frequencies than lower frequencies
12	Muscle	Power concentrated in higher frequencies than lower frequencies, but Scalp topography localizes near eyes – may be due to muscle movements during blinking.
13	Brain	Similar to 7
14	Other	Power spectrum doesn't provide strong conclusion- cannot observe noticeable peaks
15	Other	Power spectrum doesn't provide strong conclusion- cannot observe noticeable peaks and Non-dipolar scalp maps
16	Other	Similar to 14

time domain channels the spectra after removing noise components




i_LR

Normalized time domain ICA Components of the i_LR dataset






Spectrum and topography of independent components related to i_LR.



Classification of ICA components


ICA Component	Category - eye movement, muscle movement, power line noise, channel noise	Reason
1	Eye	Scalp topography strongly suggests this component captures the effects of horizontal eye movement and clearly visible large peaks in the time domain plot.
2	Eye	Similar to 1 but time domain plot doesn't show large peaks as 1.
3	Eye	Similar to 2.
4	Eye	Similar to 1
5	Eye	Similar to 2
6	Eye/Other	component is well fit by a single dipole located in the brain,
7	Brain	Scalp topography looks dipolar, Power spectrum decreases as frequency increases (1/f).
8	Brain	Scalp topography looks dipolar, Power spectrum decreases as frequency increases (1/f) most importantly a peak around 10Hz.
9	Muscle	Power concentrated in higher frequencies than lower frequencies(>20Hz).
10	Other	Power spectrum doesn't provide strong conclusion- cannot observe noticeable peaks and Non-dipolar scalp maps
11	Muscle	Power concentrated in higher frequencies than lower frequencies(>20Hz) and Scalp topography localizes near ears and neck.
12	Brain	Similar to 8
13	Other	Power spectrum doesn't provide strong conclusion- cannot observe noticeable peaks
14	Brain	Similar to 8
15	Other/Muscle	Similar to 13 but power concentrated in higher frequencies than lower frequencies(>20Hz)
16	Other	Power Spectrum consists rapid peaks, but Scalp topography looks dipolar.

time domain channels the spectra after removing noise components




i_UD

Normalized time domain ICA Components of the i_UD dataset






Spectrum and topography of independent components related to i_UD



Classification of ICA components


ICA Component	Category - eye movement, muscle movement, power line noise, channel noise	Reason
1	Eye	Scalp topographies suggest Equivalent Current Dipoles near eyes, the power concentrated at low frequencies (below 5 Hz).
2	Eye	Similar to 1
3	Other	Power spectrum doesn't provide strong conclusion- cannot observe noticeable peaks
4	Eye	Similar to 1
5	Brain	Power spectrum decreases as frequency increases (1/f) most importantly a peak around 10Hz and 20Hz.
6	Eye	Scalp topography strongly suggests this component captures the effects of horizontal eye movement and clearly visible large peaks in the time domain plot.
7	Other	Power spectrum doesn't provide strong conclusion- cannot observe noticeable peaks
8	Brain	Scalp topography looks dipolar, Power spectrum decreases as frequency increases (1/f) most importantly a peak around 10Hz.
9	Brain	Similar to 8
10	Other	Power spectrum doesn't provide strong conclusion- cannot observe noticeable peaks
11	Eye	Similar to 1 but time domain plot doesn't show large peaks as 1.
12	Brain	Scalp topography looks dipolar, Power spectrum decreases as frequency increases (1/f) most importantly a peak around 10Hz.
13	Other	Similar to 10
14	Other	Similar to 10
15	Other	Similar to 10
16	Other	Power Spectrum consists rapid peaks

time domain channels the spectra after removing noise components




teeth-blink

Normalized time domain ICA Components of the teeth-blink dataset






Spectrum and topography of independent components related to teeth-blink



Classification of ICA components


ICA Component	Category - eye movement, muscle movement, power line noise, channel noise	Reason
1	Eye	Scalp topographies suggest Equivalent Current Dipoles near eyes, the power concentrated at low frequencies (below 5 Hz) and clearly visible large peaks in the time domain plot.
2	Other	
3	Eye	Similar to 1
4	Muscle	Power concentrated in higher frequencies than lower frequencies(>20Hz) and Scalp topography localizes near ears and neck.
5	Muscle	Similar to 4
6	Muscle	Similar to 4
7	Muscle	Similar to 4
8	Eye	Similar to 1
9	Muscle	Similar to 4
10	Brain	Scalp topography looks dipolar, Power spectrum decreases as frequency increases (1/f) most importantly a peak around 10Hz, but also contains frequencies >20Hz.
11	Muscle	Similar to 4
12	Muscle	Similar to 4
13	Muscle	Similar to 4 but Scalp topography localizes near eyes.
14	Eye	Similar to 1
15	Muscle/Other	Power concentrated in higher frequencies than lower frequencies(>20Hz) and non-dipole scalp topography.
16	Other	Power spectrum consists of rapid peaks; doesn't provide strong conclusion.

time domain channels the spectra after removing noise components




noise-ON-blink

Normalized time domain ICA Components of the noise_ON_blink dataset






Spectrum and topography of independent components related to noise_ON_blink



Classification of ICA components


ICA Component	Category - eye movement, muscle movement, power line noise, channel noise	Reason
1	Eye	Scalp topographies suggest Equivalent Current Dipoles near eyes, the power concentrated at low frequencies (below 5 Hz) and clearly visible large peaks in the time domain plot.
2	Eye	Similar to 1 but time domain peaks are not visible as much.
3	Eye	Similar to 1
4	Brain	Scalp topography looks dipolar, Power spectrum decreases as frequency increases ($1/f$) most importantly a peak around 10Hz
5	Brain	Similar to 4
6	Eye	Similar to 1 but time domain peaks are not visible as much.
7	Brain	Scalp topography looks dipolar, most importantly a peak around 10Hz
8	Muscle	Power concentrated in higher frequencies than lower frequencies(>20Hz)
9	Brain	Scalp topography looks dipolar, most importantly a peak around 20Hz
10	Eye	Scalp topography strongly suggests this component captures the effects of horizontal eye movement and clearly visible large peaks in the time domain plot.
11	Brain	Similar to 4
12	Eye	Similar to 10
13	Brain	Similar to 4
14	Brain	Similar to 4
15	Brain	Similar to 4
16	Brain	Similar to 4

Time domain channels the spectra after removing noise components




noise-OFF

Normalized time domain ICA Components of the noise_OFF dataset



Spectrum and topography of independent components related to noise_OFF



Classification of ICA components

ICA Component	Category - eye movement, muscle movement, power line noise, channel noise	Reason
1	Eye	Scalp topographies suggest Equivalent Current Dipoles near eyes, the power concentrated at low frequencies (below 5 Hz) and clearly visible large peaks in the time domain plot.
2	Eye	Scalp topography strongly suggests this component captures the effects of horizontal eye movement and clearly visible large peaks in the time domain plot
3	Eye/ Brain	Similar to 1 but a peak around 10Hz
4	Other/Channel noise	Very focal scalp topography.
5	Other	Power spectrum doesn't provide strong conclusion, non-dipole scalp topography.
6	Brain	Scalp topography looks dipolar, Power spectrum decreases as frequency increases ($1/f$) most importantly a peak around 10Hz
7	Brain	Similar to 6
8	Brain	Scalp topography looks dipolar most importantly a peak around 10Hz
9	Other	Similar to 5
10	Brain	Similar to 6
11	Other	Similar to 5
12	Brain	Similar to 6
13	Other	Similar to 5
14	Other/Brain	Peak around 10 Hz
15	Brain	Power spectrum decreases as frequency increases ($1/f$), A large peak around 10 Hz
16	Other	Power spectrum consists of rapid peaks; doesn't provide strong conclusion.

time domain channels the spectra after removing noise components

